

myFUGETM Mini

Mini Centrifuge



Eight Tube Capacity

The myFugeTM Mini centrifuge is a personal centrifuge that virtually fits in the palm of your hand. Ideal for quick spin downs of microtubes and PCR tubes, the MyFuge is extremely easy to use. Simply close the lid and the rotor instantly reaches 6,000 rpm. Open the lid, and the improved braking system quickly and smoothly decelerates the rotor to a complete stop in just one second.

An eight-position microtube rotor is included and comes installed. Also included is a second rotor with increased capacity for four 0.2ml PCR strips (or 32 x 0.2ml tubes) which is stored in a convenient compartment on the bottom of the centrifuge. This compartment is also used for securing the included tube adapters for 0.5ml tubes.

A new **Clik-n-LokTM** rotor attachment system eliminates the requirement for set screws and tools during rotor exchange. Simply remove and attach rotors by hand, pulling or pressing the rotor onto the motor shaft.

- + **Clik-n-LokTM** Instant rotor attachment
- + **Eight position microtube rotor**
- + **Storage compartment for PCR rotor and 0.5ml & 0.2ml adapters**
- + **Decelerates in just 1 second**
- + **Four position PCR strip rotor**

Technical Data

Speed:	6,000rpm / 2,000 xg
Capacity:	8 x 1.5 / 2.0ml tubes 32 x 0.2ml PCR tubes 4 x PCR strips (8x0.2ml)
Dimensions:	5.7 x 6.9 x 4.5 in 14.5 x 17.8 x 11.4 cm
Weight:	2.2 lbs / 1 kg
Electrical:	100 - 240V, 50-60Hz



Storage Compartment



Includes Both Rotors

Ordering Information

C1008-B	MyFUGE Mini Centrifuge, blue lid
C1008-C	MyFUGE Mini Centrifuge, clear lid
C1008-G	MyFUGE Mini Centrifuge, neon green lid
C1008-P	MyFUGE Mini Centrifuge, purple lid
C1008-R	MyFUGE Mini Centrifuge, red lid

*All include microtube rotor, PCR rotor, 0.5ml & 0.2ml adapters

Benchmark
Scientific

Benchmark Scientific Inc.
PO Box 709
Edison, NJ 08817
Phone: 908-769-5555
Fax: 908-222-1864

Web: www.BenchmarkScientific.com
Email: Info@BenchmarkScientific.com