




EXTRA SECURITY

FLEXIBLE STORAGE

OPTIONAL ACCESSORIES

MAIN DIMENSION FOR PRE 120





UNIT INFORMATION

Item code	DAI 6040
Dimensions exterior WxDxH (mm)	600x600x835
Controller (name)	Carel
Capacity (L)	117
Dimensions interior WxDxH (mm)	480x410/240x690
Insulation (mm)	50
Temperature range (°C)	+2/ +8
Max. ambient temperature (°C)	32
Cooling technology	Air forced cooling

POWER

Power supply (V)	230
Frequency (Hz)	50/60
Power consumption (KWH/24H)	1,2
AMP (A)	0,5
Power (Watt)	107
Horse power (HP)	1/6
Battery backup for alarms and temperature display in the event of power cut (approx hours)	no

ALARM FUNCTIONS

Alarm display as text - not codes	no
Visual / Acoustic alarm	yes
Power failure alarm	no
Adjustable high/low temperature alarm	yes
Open door alarm	yes
Probe failure alarm	yes
Low battery alarm	no
Compressor alarm	no

LOGGER & EXT. CONNECTION

Contact for remote alarm	yes
Prepared for GSM alarm module	yes
Data logger	no
RS 232 interface	no
USB for software upload	no
USB readout for logged data	no
Chart recorder	no

DISPLAY FEATURES

Digital display	yes
Defrost manual / auto	yes
2 extra probes can be mounted inside	no
Password protection for turning unit on/off	no
Temperature graph on display	no

OTHER EQUIPMENT

Porthole for external probe etc.	yes
Porthole position	back
Porthole diameter	11
Temperature sensor	NTC
Lighting	yes
Key lock / Card lock	yes / no
Castors / with break in front	2 small back / no
Shelves (pcs.) / drawers (pcs.)	2+1 / no
Adjustable shelves (pcs.) / drawers (pcs.)	3 / no
Max load per shelf	2 - 15kg & 1 - 10
Vaccum release port	1
Basket (pcs.)	no
Sub lids / sub doors (pcs.)	yes
Reversible doors (pcs.)	yes

REFRIGERATION & CABINET

No. of compressors	1
Refrigerant	R600a
Inner liner material	Coated steel
Cabinet material	Ccoated steel
Colour (exterior)	white
Foam agent	polyurethane

SHIPPING

Weight (KG)	51
Shipping weight (KG)	59
Shipping dimensions WxDxH (cm)	61,2x80x99,5
Shipping volume (m³)	0,49

Arctiko A/S

Lammefjordsvej 5 · DK-6715 Esbjerg N · Phone: +45 7020 0328
info@arctiko.com ■ arctiko.com

