

AmpliPur™

EXPERT TIPS



EXCELLENT FIT FILTER TIPS FOR MOLECULAR BIOLOGY APPLICATIONS

| Excellent Fit

| Ideal for Molecular Biology Applications

| ISO 8655 Compliant



 **GILSON®**



Excellent Fit*

AmpliPur Expert Tips are designed to fit various pipette brands for maximum convenience on the bench. Equipped with a filter, they also help prevent cross-contamination.


Ideal for Molecular Biology Applications

These tips significantly contribute to molecular testing performance by ensuring sample integrity during the pipetting process.

ISO 8655 Compliant

Manufactured according to the certified Gilson quality standard, AmpliPur Expert Tips comply with ISO 8655 requirements and are certified free of RNase and DNase, preserving your sample integrity.

*Pipettes compatibility need to be verified by the user.



Technical Information	
Gilson Quality System	ISO 9001
Product Compliance	ISO 8655
Analysis Certificate	Free of RNase and DNase Hazardous Materials (2015/863/EU)
Materials	Tips: Polypropylene (PP)
	Filter: Polyethylene (PE)
	Rack: Polypropylene (PP)
	Heat shrink: Plastic Optical Fibre (POF)
Packaging	10 racks of 96 tips (960 tips)
Disposable	Tips: Yes
	Rack: Reusable
Warranty	3 years

AmpliPur™ Expert Tips		
Volume Range	Nominal Volume of Compatible Pipettes**	Part Number (Commercial Units)
0.5-10 µL	2 µL / 10 µL	F174101
1-20 µL	20 µL	F174201
10-200 µL	100 µL*** / 200 µL 300 µL	F174301
100-1000 µL	1000 µL	F174401
100-1200 µL	1000 µL / 1200 µL	F174501

**Within ISO 8655 specifications.

***Compatible with many of 100 µL pipettes, except Gilson P100 single channel.

AmpliPur Expert Tips are only available in select countries. Please contact your local Gilson representative for more information.



gilson.com/contactus

ISO 9001
CERTIFIED

LT800607/A | ©2021 Gilson S.A.S.
All rights reserved. All trademarks
are the property of Gilson, Inc.
and its subsidiaries.





Click on the symbols     to go directly to the product pages.

Purity Grade	Pipette Compatibility	Filter	Packaging Type (Code Description)	Volume Range															
				0.1-10 µL	0.1-20 µL	1-20 µL	2-30 µL	1-100 µL	10-100 µL	1-200 µL	2-200 µL	20-200 µL	1-300 µL	10-300 µL	100-1000 µL	100-1200 µL	50-1250 µL	0.5-5 mL	1-10 mL
PIPETMAN DIAMOND TIPS Sterile SAL 10⁻⁶* - Free of detectable DNA, RNA, DNase, RNase, ATP, endotoxins (pyrogens) - Trace metal release tested	PIPETMAN	With filter	TIPACK (DFxxxST)	●			●■			●■			●■□		●■	●			
			BLISTER REFILL (DFxxxST)	●◆			●◆			●◆			●◆	●◆	●◆				
			BLISTER REFILL - Certified 384-well (DSFxxxST)	◆			◆				◆■□			◆					
		Without filter**	TIPACK (DxxxST)	●	●						●■□			●■	●	●			
			TOWERPACK (DxxxST)	●	●						●			●	●				
			STERILPACK (DxxxST)								●				●				
			BLISTER REFILL (DxxxST)		●◆						●◆			●◆	●◆	●◆			
BLISTER REFILL - Certified 384-well (DSxxxST)		◆						◆■□				◆							
PIPETMAN EXPERT TIPS Sterilized SAL 10-3* Free of human and bacterial DNA, DNase, RNase, non pyrogenic, no inhibition of real time PCR <i>Available in select countries only</i>	Most leading brands	With filter	TIPACK (EFxxxST)	●		●		●				●			●				
		Without filter**	TIPACK (ExxxxST)	●	●					●			●			●			
			TIPACK - Gel loading flat or round tips (ExxxxST)**							●									
			MEGAPACK (ExxxxST)		●					●									
PIPETMAN DIAMOND TIPS Autoclavable** - Manufactured in clean room environment - Trace metal release tested	PIPETMAN	Without filter	TOWERPACK (Dxxx)	●	●					●			●	●					
			TIPACK (Dxxx)	●	●					●■□			●■	●	●		●	●	
			BLISTER REFILL (Dxxx)		●◆					●◆			●◆	●◆	●◆				
			BLISTER REFILL - Certified 384-well (DSxxx)		◆					◆■□				◆					
PIPETMAN DIAMOND TIPS Autoclavable** Trace metal release tested	PIPETMAN	Without filter	EASYPACK (Dxxx)	●	●					●■□			●■	●	●		●	●	
			RELOAD PACK (Dxxx)		●					●■□				●					
			ECOPACK (Dxxx)	●	●					●■□			●■	●	●				
PIPETMAN EXPERT TIPS Autoclavable** Standard <i>Available in select countries only</i>	Most leading brands	Without filter	TOWERPACK (Exxx)**		●				●			●			●				
			TIPACK (Exxx)	●	●				●			●			●				
			TIPACK - Gel loading flat or round tips (Exxx)**	●					●										
			EASYPACK (Exxx)		●				●										

PIPETMAN® DIAMOND Tips are designed to fit PIPETMAN® perfectly and guarantee maximum accuracy and precision in conjunction with your Gilson pipette.

PIPETMAN® EXPERT Tips can be used with most leading pipette brands.

*Sterility Assurance Level SAL (ISO 11137): quantitative value to assure sterility. SAL 10⁻⁶: probability of 1 in 1.000.000 of finding a non-sterile unit.

**Autoclavable: tips are fully autoclavable except gel loading tips. Autoclaving does not destroy traces of DNA and RNases. Sterile not filtered PIPETMAN DIAMOND Tips and Sterilized not filtered PIPETMAN EXPERT tips are autoclavable.

TOWERPACK Tips are autoclavable, but the plastic wrap and the tray of the TOWERPACK are not. We recommend to place the racks in the tip boxes before autoclaving.

Empty reload boxes and racks are available for RELOAD, BLISTER REFILL, TOWERPACK, ECOPACK, EASYPACK packaging options.